



In-booth access to audiometric tools

Skip the external monitor



Screen Mirror lets audiologists view and interact with audiometric data directly inside the booth.

- Review completed audiograms with patients inside the booth.
- View a video otoscope directly on a Flex VRA monitor.
- Support counseling with on-screen annotation, drawings, and visuals.
- Control audiometric software during testing with a wireless keyboard.

Reduce room exits, improve counseling flow, and keep the patient engaged.

Designed and assembled in Austin, Texas.

Signal Transmission	HDMI
Latency	Instantaneous
Connection	HDMI Video Caputure to USB
HDMI Length	15 foot
Warranty Period	2 Years
Flex SW Requirements	2.0 or higher

Combination keyboard and trackpad included for play audiometry in the booth.

