



IN-BOOTH ACCESS TO AUDIOMETRIC TOOLS

Skip the external monitor.

Screen Mirror lets audiologists view and interact with audiometric data directly inside the booth, on a Flex VRA monitor.



- ✓ Review completed audiograms with patients inside the booth
- ✓ Support counseling with on-screen annotation, drawings and visuals
- ✓ View a video otoscope directly on a Flex VRA monitor
- ✓ Control audiometric software during testing with a wireless keyboard

Reduce room exits, improve counseling flow, and keep the patient engaged.

Specifications

Video input	HDMI
Latency	Instantaneous
Connection	HDMI video capture to USB
HDMI cable	15 ft (included)
Software	Flex VRA software 2.0 or higher
Warranty	1 year

FLEX-SCREENMIRROR



Includes the Flex video input, a 15 ft HDMI cable and a combination keyboard and trackpad for use in the booth.