


# FLEX TOY



## TRADITIONAL VRA

Available as a stand-alone system or integrated into Flex VRA video systems

### WHY FLEX TOY

- Backlit animated figures with sound that include over 15 different nursery rhythms
- Startle-free silent spinning motor (Big Box)
- Distract Mode (activate to engage patient during Tymp/OAE/Otoscopy)
- Can be controlled by computerized audiometer or GSI Audiostar Pro
- Adjustable light color
- Customizable toy and sound options



Custom rubber coated toys designed to be inviting, novel, and interesting for pediatric patients

### MORE INFORMATION

512-777-0016

sales@flexvra.com

www.flexvra.com/flextoy

www.youtube.com/flexvra

www.fb.com/flexvra

@flexvra


### OPTIONAL ADD ONS



FlexPC controls Flex Toy using a computerized audiometer



Control Flex Toy directly on the GSI Audiostar Pro



Integrate Flex Toy into any Flex Video system to access up to 6 reinforcers and simultaneously activate toy and video reinforcers

# STAND-ALONE CONFIGURATIONS



## FLEX-TOY-BBB

Left, right, and center big boxes with lights, motion and sound.



## FLEX-TOY-BsB

Left and right big boxes with lights, motion and sound. Small center with no motion.



## FLEX-TOY-sBs

Left and right small box with lights, sound, and no motion. Big center box with motion/sound



## FLEX-TOY-sss

Left, right and center small boxes with lights, sound and no motion.



Give Flex Toy a try at [flexvra.com/demo](http://flexvra.com/demo)

## THE REMOTE



## SPECIFICATIONS

**Remote Signal:**  
Infrared

**Power Req:** 12 V

**Big Box Dimensions:**  
10x10x6in

**Small Box Dimensions.:**  
6x5.5x4in

**Warranty.:** 1 year warranty from date of shipment. 60 day free return from date of shipment.

