



## Visual Reinforcement Operant Conditioning Audiometry



## Better engagement. Better Results.

The Flex VROCA button enables patient-initiated responses that trigger visual reinforcement. The wireless and battery-free design allows testing anytime and anywhere without added setup, and the built-in compatibility with Flex VRA gives clinicians precise control over the reinforcer. Eliminating cables reduces distraction and allows children to respond naturally without shifting posture or attention.

**The result is improved response consistency and greater confidence in behavioral test results.**

### Specifications

Designed for fast, reliable clinical use

Signal Transmission	Radio Frequency
Activation Latency	.5 ms (Instantaneous)
Button Power	Kinetic (Button Press Sends the Charge)
USB Receiver Power	USB 5V
RF Frequency	915 MHz
Weight	.25 lbs
Button Dimensions	4" (L) x 4" (W) x 4" (H)
Warranty Period	2 Years

